
yuwell

DE

GB

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Thank you for your trust in our product!

TABLE OF CONTENTS

For safe use of your new product, please pay particular attention to the following instructions for use and safety information.

If you have any questions, please contact our customer service team. You can reach our service staff on the toll-free hotline: +49 (0)6421 - 30 30 800

This user manual has been compiled based on the best information currently available. The information in this document is subject to change. The product you have purchased may differ slightly from the product described in the user manual. The latest user manual can be requested from the online shop. It is regularly updated to provide you with the latest information.

- Before using the device for the first time, be sure to read these instructions carefully and follow the safety instructions to ensure proper use.
- Any liability is excluded in the event of failure to observe the safety rules specified in the operating instructions.
- Keep these instructions for future reference and pass them on to other users.
- Failure to follow these operating instructions and/or observe the safety measures may result in falls and serious injuries.

 **NOTE:** This English translation is based on the German version.

1. IMPORTANTWICHTIGE HINWEISET NOTES	page 42
1.1 Explanation / Symbols in the instructions for use	
1.2 Explanation of type plate	
1.3 Barrier-free instructions for use	
1.4 Importer	
2. INSTRUCTIONS FOR USE	43
2.1 Declaration of conformity	
2.2 Intended purpose	
2.3 Indications	
2.4 Contraindications	
2.5 Patient target group	
2.6 Target group	
2.7 Safety instructions	

3. FUNCTIONS	page 52
3.1 Intended use	
3.2 Technical safety features	
3.3 Specifications	
4. HANDLING	53
4.1 Unpacking	
4.2 Inspection	
4.3 Storage	
5. OPERATION AND INSTALLATION	54
5.1 Functions	
5.2 Getting started	
5.3 Oxygen-absorbing operation	
5.4 Alarm signals	
5.5 Timer function	
5.6 Symbols	
5.7 Adjusting the brightness	
5.8 Switch off	
5.9 Accessories	
6. MAINTENANCE AND CARE	65
6.1 Cleaning	
6.2 Clean or replace filters	
6.3 Clean humidifiers	
6.4 Cleaning the fire safety valve	
6.5 Replace fuse tube	
6.6 Check the alarm system	
6.7 Instructions for preparation and reprocessing	
7. TROUBLESHOOTING	69
8. FURTHER IMPORTANT INFORMATION	72
8.1 Functional diagram of oxygen production	
8.2 Electrical function diagram	
8.3 Package contents	
8.4 Disposal of equipment	
8.5 Language	
9. EMC INFORMATION	74

GB


1. IMPORTANT NOTES

1.1 Explanation / Symbols in the instructions for use

-  Follow the instructions for use!  Manufacturer
-  Order number  Date of manufacture
-  Batch number  European authorised representative
-  Importer  Swiss authorised representative
-  The instructions for use contain important safety-related information
-  This product is in conformity with Regulation (EU) 2017/745
-  This product is a Class IIb medical device
-  **Allergy warning:** The high-quality rubber handles and wheels contain natural rubber latex. Natural rubber latex can trigger allergic reactions, including anaphylactic shock.
-  Do not use the product at temperatures below -35 °C or above +60 °C
-  This product does NOT belong in household waste
-  This product is designed for indoor and outdoor use

1.2 Explanation of type plate

- | | |
|--------------------------------------|---------------------------------------|
| 1. Manufacturer | 5. SN (serial number) |
| 2. Model | 6. manufacturer's date |
| 3. UDI-DI (Unique Device Identifier) | 7. symbols (see chapter 1.1 and 5.) |
| 4. total weight | 8. European authorised representative |

1.3 Barrier-free instructions for use

-  These instructions for use contain all the necessary information for adjusting and operating the product. The latest edition of the instructions for use is available online as a PDF file or can be requested from our customer service department. Enlarged formats can be provided for visually impaired persons. Blind patients must be instructed by their accompanying persons.

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2. INSTRUCTIONS FOR USE

2.1 Declaration of conformity

This oxygen concentrator is intended for oxygen supplementation. The oxygen concentrator is a Class IIb medical device and complies with the requirements of the European Medical Devices Regulation (EU) 2017/745. A corresponding EU declaration of conformity is available.

2.2 Intended purpose

This oxygen concentrator is intended for oxygen supplementation.

2.3 Indications

For the administration of oxygen as part of long-term oxygen therapy for hypoxaemia requiring treatment – at least 15 hours per day – in a home environment.

2.4 Contraindications

Users with oxygen toxicity or oxygen allergy should NOT use this oxygen concentrator. This device is intended as supplemental oxygen and is NOT intended as a life-sustaining device. Users who require continuous oxygen must have alternative backup sources of power and oxygen in case of power failure or loss of oxygen supply.

2.5 Patient target group

The product may only be used in adult patients and is not suitable for children or adolescents.

2.6 Target group

-  **WARNING:** The product is intended for use by both medical professionals and lay persons. The patient himself is considered the intended user of the product.

2.7 Safety instructions

Symbol	Description
WARNING 	Describes a hazard or unsafe practice which, if not avoided, could result in serious injury, death or property damage.
CAUTION 	Describes a hazard or unsafe practice which, if not avoided, could result in minor injury or property damage.

SAFETY INSTRUCTIONS

I. Important notice

- To ensure your safety and the proper functioning of the device, please observe the following instructions before using the oxygen concentrator:
- Risk of fatal electric shock: Never open or disassemble the oxygen concentrator yourself. Repairs may only be carried out by qualified service personnel.
- Before first use: Read this user manual carefully and in its entirety before using the device.

II. Before installation

- The oxygen concentrator should always be kept upright to prevent damage during transport.
- If the mains voltage is unstable and outside the normal voltage range, please use the voltage stabiliser.
- Use a suitable, safe power supply unit and a suitable connection box.

SAFETY INSTRUCTIONS

- The housing of the oxygen concentrator must not be dismantled by untrained persons.
- Removing the oxygen concentrator or replacing internal components may result in injury to personnel or damage to the device.

III. Place of use

- You are free to decide in which room of your home you would like to use the oxygen concentrator. Thanks to its castors, the oxygen concentrator can be easily moved from one room to another.
- Ensure that the device is at least 10 cm away from walls, curtains, furniture or other obstacles on all sides.
- Ensure that the air inlet and outlet of the oxygen concentrator are in a well-ventilated environment.
- The oxygen concentrator should be positioned so that it is protected from pollutants or vapours.
- The oxygen concentrator must be protected from heat sources, fire sources, moisture and environments with extremely high or low temperatures.
- No objects or containers may be placed on the oxygen concentrator.
- Always position the oxygen concentrator so that the operator can hear the audible alarm.
- Do not use the oxygen concentrator in an MR (magnetic resonance) environment.

IV. Use

- Warning: To ensure that you receive the therapeutic amount of oxygen appropriate for your health condition, the oxygen concentrator must not be used without supervision.
 - Only be used with the settings that have been individually determined for you or prescribed for your activity level with accessories.
 - Only be used with the specific combination of parts and accessories that meet the specifications of the concentrator or accessory manufacturer.
- Warning: Using the device at altitudes above 2,000 metres, at temperatures below 5 °C or above 32 °C, or at humidity levels above 90% may affect the oxygen level and airflow. This may reduce the effectiveness of the therapy.
- The oxygen supply settings on the oxygen concentrator should be checked regularly to ensure the effectiveness of the therapy.

GB

SAFETY INSTRUCTIONS

- The oxygen supply setting must be determined individually for each patient, taking into account the configuration of the device to be used, including accessories.
- Only lubricants recommended by the manufacturer may be used.
- No parts, accessories or adapters that are not approved by the manufacturer may be used.
- The oxygen concentrator must not be used in parallel or in series with other oxygen concentrators or oxygen therapy devices.
- Avoid using the device while bathing. If continuous use is required according to a doctor's prescription, please set up the oxygen concentrator in another room, at least 2.5 metres away from the bathroom.
- The device may only be used in accordance with a doctor's prescription and the instructions in this manual.
- If the patient or a caregiver feels at any time that the oxygen supply is insufficient, the attending physician and/or the manufacturer must be contacted immediately.
- The flow rate may only be adjusted on the express instruction of a doctor.
- It takes 15 minutes for the oxygen concentrator to start up and reliably deliver the set flow rate and oxygen concentration.
- To achieve optimum performance, the oxygen concentrator should not be opened or switched off frequently, but must be restarted after 3 to 5 minutes. Shorter operating times can shorten the maximum service life of the product.
- The service life at 8 hours per day is 3 years.

V. Maintenance

- The oxygen concentrator has been designed to require routine preventive maintenance every 2 years. This work must only be carried out by qualified personnel, including medical professionals or authorised, factory-trained service technicians who are familiar with the maintenance and performance adjustment of the device.
- Warning: Maintenance and servicing work must not be carried out while the oxygen concentrator is in operation.
- The manufacturer will provide authorised service personnel with circuit diagrams, parts lists, calibration instructions and other technical information upon request. These documents are intended for the professional repair of components that have been classified as repairable by the manufacturer. Please leave all repair and maintenance work to qualified personnel only to ensure the safety and functionality of the device.

SAFETY INSTRUCTIONS

VI. Funkstörung

- This device has been tested and complies with the electromagnetic compatibility limits specified in IEC 60601-1-2. These limits are intended to provide reasonable protection against electromagnetic interference in a typical medical facility.
- In other facilities, even the low levels of electromagnetic emissions permitted by the above standards may cause interference. To determine whether the emissions from the concentrator are causing the interference, turn off the concentrator. If the interference from the other device(s) stops, the concentrator is the source of the interference. In these rare cases, the interference can be reduced or eliminated by one of the following measures:
 - Position, move or increase the distance between the devices.
 - Connect the device to a power outlet that belongs to a different circuit than the other device(s).
 - Contact the manufacturer or maintenance personnel.

VII. Avoid risks of burns, electric shocks, fire and injuries

- Oxygen enrichment poses a fire hazard. Do not use the oxygen concentrator and accessories near sparks or open flames.
- Before and during oxygen therapy, use only oxygen-compatible water-based lotions or ointments.
- Never use petroleum- or oil-based lotions or ointments to avoid the risk of fires and burns.
- Do not lubricate any connections, joints, hoses or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
- Use only manufacturer-recommended replacement parts to ensure proper operation and avoid the risk of fire and burns.
- Oxygen promotes the formation and spread of fires. Do not leave the nasal cannula or mask on bed covers or chair cushions when the oxygen concentrator is turned on but not in use. The oxygen makes these materials more flammable. Turn off the oxygen concentrator when not in use to prevent oxygen enrichment.
- Smoking during oxygen therapy is dangerous and can cause burns to the face or death. Smoking and open flames are prohibited in the room where the oxygen concentrator or oxygen-carrying accessories are located. When smoking, always turn off the oxygen concentrator, remove the cannula and leave the room where the cannula, mask or oxygen concentrator is located. If you cannot leave the room, you must wait 10 minutes after turning off the oxygen concentrator.

SAFETY INSTRUCTIONS

- Open flames are dangerous during oxygen therapy and can cause fires or death. Open flames must not be brought closer than 2 metres to the oxygen concentrator or oxygen-carrying accessories.
- During oxygen therapy, special care must be taken to reduce the risk of fire. Smoking is not permitted while using the device. Matches, lit cigarettes or other sources of ignition must not be present in the room where the device is located. No smoking signs must be clearly displayed. Textiles and other materials that do not normally burn will ignite easily and burn intensely in oxygen-enriched air. Failure to heed this warning may result in serious fire, property damage, personal injury, or death.
- When oils, fats, or greasy substances come into contact with pressurised oxygen, violent spontaneous combustion may occur. These substances must be kept away from the oxygen concentrator, hoses and connectors, and all other oxygen equipment.
- Avoid sparking near the oxygen concentrator. This also applies to sparks caused by static electricity or friction.
- If you feel unwell or experience a medical emergency during oxygen therapy, seek medical attention immediately to prevent harm.
- Geriatric, paediatric or other patients who are unable to communicate their discomfort may require additional monitoring to relay information about the discomfort and/or medical emergency to the appropriate nursing staff to prevent harm.
- Careful supervision is required when the oxygen concentrator is used near children or physically disabled persons.
- Please note that the power cord and nasal cannula may pose a tripping or strangulation hazard. Always lay the power cord and nasal cannula so that they cannot be trapped by castors or other objects.
- Ensure that a child does not choke on a small part that has come loose from the oxygen concentrator.
- Do not touch the oxygen concentrator with wet hands.
- The oxygen concentrator must not be set up or stored in such a way that it can fall into water or other liquids.
- Do not reach for an oxygen concentrator that has fallen into water. Disconnect the device from the mains immediately.
- Never leave the oxygen concentrator unattended when it is connected to a power outlet.
- Do not make any modifications to the oxygen concentrator. Modifications may result in hazards for the user.

SAFETY INSTRUCTIONS

- This oxygen concentrator and its accessories do not contain natural rubber latex.
- The oxygen concentrator and accessories do not contain phthalates.
- Prolonged contact with application parts or other accessories does not cause skin irritation.
- Do not move the oxygen concentrator when it is switched on.
- Oxygen therapy can be dangerous in certain circumstances. The manufacturer recommends consulting a doctor before use.
- To reduce the risk of infection to other users or operators when reusing the oxygen concentrator, the housing should be cleaned with a mild household cleaner before reuse. All external tubing, cannulas, humidifiers or other accessories should be replaced.
- The dealer or responsible person should instruct the operator to determine the patient's need for supplemental oxygen in the event of a failure of the oxygen concentrator or power supply:
 - a) during installation, based on
 - the patient's condition,
 - the environment in which the patient lives,
 - the possibility of supplying the patient with additional oxygen, and
 - b) at regular intervals, as these characteristics change.Lay users should contact the manufacturer or the manufacturer's representative:
 - if necessary, to obtain assistance with the installation, operation or maintenance of the device, or
 - to report unexpected operating conditions or incidents.
- The oxygen concentrator's over-temperature alarm can help limit the spread of fire in the event of ignition.
- All parts of this oxygen concentrator are intended for use in the patient environment.
- Lint, dust, animal hair and pests can clog the oxygen concentrator's air inlet and outlet. Check and clean it regularly.
- The heat from a fireplace or radiant heater can accelerate the ageing of the oxygen concentrator's internal components.
- Moisture from a nebuliser or steam boiler can accelerate the ageing of the molecular sieves in the oxygen concentrator.
- Ensure that children do not change the settings on the oxygen concentrator.

3. FUNCTIONS

This user manual contains information about your concentrator and serves as a reference guide for using the concentrator. Please read it carefully and keep it in a safe place.

3.1 Intended use

- **Intended medical use:** This oxygen concentrator is intended for use as an oxygen supplement device in professional healthcare facilities and in the home. It delivers a high concentration of oxygen to individuals who require oxygen therapy.
- **Intended patient group:** Adults only.
- **Intended body part or type of tissue with which it comes into contact or interacts:** During oxygen inhalation, the nasal cannula is in direct contact with the user's facial skin and nasal cavity.
- **Intended user profile:** Healthcare professionals or lay persons; the patient is also the intended user.
- **Usage environment:** Use in hospital or at home.
- **How it works:** The oxygen concentrator uses air as a medium and a molecular sieve as an adsorbent to produce oxygen through pressure swing adsorption.

3.2 Technical safety features

The oxygen concentrator has been equipped with a variety of well-designed features to ensure safe, reliable and comfortable use:

- Robust plastic housing – for maximum safety and durability.
- Time accumulation function – displays the total usage time directly on the screen.
- Time measurement can be switched off – allows the time measurement to be reset or paused if necessary.
- Integrated pressure relief valve on the compressor – provides additional safety during use.
- Alarm function in case of power failure and loss.
- Alarm function in case of device failure (including pressure failure, compressor failure, low oxygen concentration, low oxygen flow, overtemperature).
- Compressor with overheating protection – protects both the compressor and the oxygen concentrator and ensures long-term operational safety.

3.3 Specifications

- Power supply: 230 V, 50 Hz
- Input power: 350 VA
- Power consumption in standby mode: < 2 W
- Recommended maximum flow rate: 5 l/min
- Maximum power consumption: 360 W
- Oxygen concentration at a nominal outlet pressure of 0 kPa (measured after a 15-minute warm-up period): 5 l/min: 94.4%
- Maximum permissible pressure: 70 kPa
- Sound pressure level (measured at a distance of 1 m from the device) 40 dB(A) and sound power level 53 dB(A) at 3 l/min and 5 l/min with an uncertainty of 2 dB(A), measured according to the noise measurement method in accordance with ISO 80601-2-69 using the basic standard ISO 3744. 7.

3. FUNCTIONS

- Audible alarm: 48 dB(A) or more for failure alarm; 40 dB(A) or more for power failure.
- The nominal range of both oxygen input flow and oxygen concentration as a function of flow:
 - Tested under STPD conditions (101.3 kPa, 20 °C, dry) and under the nominal operating conditions specified in the manual. (Fig. 1)

Flow rate	Oxygen Concentration
1L/min	94,1 %
2L/min	94,1 %
3L/min	94,0 %
4L/min	93,9 %
5L/min	94,4 %

Fig. 1

- The oxygen concentration may be affected by exceeding the nominal ranges of ambient temperature, humidity and air pressure. The measurement uncertainty of the flow rate is $\pm 10\%$. The measurement uncertainty of the oxygen concentration is $\pm 3\%$.
- Altitude: no more than 2000 m above sea level.
- Net weight: 17.5 kg
- Dimensions: 40×30×55 (cm)
- Working system: continuous
- Minimum operating time: 15 minutes
- Electrical classification: Class II device, type BF application part, IP21
Application part: Nasal cannula
- Electrical classification: Overvoltage category: II; Pollution degree: 2; Altitude: ≤ 2000 m
- Safety system:
 - Power failure: Alarm and shutdown
 - Pressure failure: Alarm and shutdown
 - Compressor failure: Alarm and shutdown
 - Low oxygen concentration: Alarm
 - Low flow alarm: Alarm and shutdown
 - Overtemperature alarm: Alarm and shutdown
- Normal operating conditions (with oxygen concentration status display)
 - Temperature range: 5 °C to 32 °C
 - Relative humidity: 15% to 90%, non-condensing
 - Atmospheric pressure: 86 kPa to 106 kPa

GB

3. FUNCTIONS

 **WARNING:** Operating conditions outside the rated ranges of ambient temperature, humidity and air pressure may impair the performance of the oxygen.

- Oxygen outlet temperature: $\leq 46\text{ °C}$
- Temperature of the application part (nasal cannula): $\leq 41\text{ °C}$
- Cannula length: NOT more than 15.2 m and not twisted.
- Storage and transport conditions:
 - Temperature range: $-20\text{ °C} \sim 60\text{ °C}$
 - Relative humidity: $\leq 93\%$, non-condensing.

 • **CAUTION:** The device should be stored away from strong sunlight and corrosive gases, in a well-ventilated indoor area. The device may only be transported and used in a vertical position.

 • **CAUTION:** If the oxygen concentrator has been stored at very high or very low temperatures, it must be left at room temperature (approx. 20 °C) for at least 4 hours before use. Only then is it ready for operation and can be used safely.

4. HANDLING

4.1 Unpacking

 **Caution:** Containers and packaging materials should be kept until the oxygen concentrator is used.

1. Check the box or other packaging for obvious damage. If damage is apparent, please notify the carrier or your local dealer.
2. Remove all loose packaging materials from the box.
3. Carefully remove all components from the box.

4.2 Inspection

1. Inspect the exterior of the oxygen concentrator for nicks, dents, scratches or other damage.
2. Check all components.

4.3 Storage

1. Store the repackaged oxygen concentrator in a dry place.
2. Do not place any other objects on top of the oxygen concentrator.

GB

5. OPERATION AND INSTALLATION

5.1 Functions

Fig. 2: Front view




Fig. 3: Rear view


5. OPERATION AND INSTALLATION

Fig. 4: Control panel

GB

Fig. 5: Components of the humidifier


5. OPERATION AND INSTALLATION

5.2 Getting started

NOTE: Before use, inspect the power cord/plug and the exterior of the oxygen concentrator for nicks, dents, scratches or other damage. If necessary, have the device inspected and repaired by qualified service personnel.

1. Unscrew the bottle from the humidifier by turning it counterclockwise. Fill the bottle with clean water (or distilled water) up to the mark between MAXIMUM and MINIMUM. Do not fill the humidifier bottle above the MAXIMUM mark. (Figure 5)
2. Screw the bottle clockwise until tight. (We recommend using the Yuwell humidifier, whose preferred position is shown in Figure 2)
3. Connect the power supply.
4. If you want to move the oxygen concentrator, unlock the locks on the four castors. (Figure 6)

Fig. 6: Castors

- CAUTION:**
1. The power cord of the oxygen concentrator is not detachable. If the power cord is damaged, contact the service personnel to have it replaced..
 2. Keep the power cord away from heated surfaces.
 3. Do not move or relocate the oxygen concentrator by pulling on the power cord.
 4. Do not use extension cords with this device.
- NOTE:** The oxygen concentrator can be used during the initial warm-up period (approximately 15 minutes) while the oxygen concentration rises to its maximum value.

5. OPERATION AND INSTALLATION

5.3 Oxygen-absorbing operation

Switch on:


When the power button is pressed to the 'I' position, the message 'HELLO' appears on the display and the blue, green and yellow indicators light up simultaneously, indicating that the oxygen concentrator is functioning properly. After 2 seconds, the oxygen concentrator switches to standby mode and only the green indicator remains lit. Press the standby button, the display backlight will turn on and show the time and total operating hours, then the oxygen concentrator will switch to normal operating mode. When the oxygen concentrator is in operation, it will emit 'clicking' sounds every few seconds, which is the normal reverse and exhaust sound.

Flow rate:

Set the flow meter knob to the desired flow rate (the measurements should be based on the black centre of the float). Turn the flow meter knob upwards to increase the flow rate and downwards to decrease the flow rate (Figure 7). In the meantime, the humidifier should have air bubbles around the filter element. Then oxygen comes out of the oxygen outlet.


Connect the fire safety valve to the humidifier outlet in the direction and position shown in Figure 8. Connect the nasal cannula to the fire safety valve, and the other end is worn by the patient, and oxygen inhalation can begin.

Fig. 7: Set flow rate

5. OPERATION AND INSTALLATION

Fig. 8: Fire safety valve

 **WARNING:** For the duration of oxygen inhalation and adjustment of oxygen flow, please follow the doctor's instructions.

NOTE: If the flow rate on the flow meter ever falls below 0.5 l/min, check whether the tubes or accessories are blocked or kinked, or whether the humidifier is defective.

NOTE: Connect the nasal cannula to the oxygen outlet port of the oxygen concentrator. With the oxygen concentrator turned on, set the flow meter to the desired flow rate. Air should flow freely to the nasal cannula. You should be able to hear or feel the oxygen flow to the prongs of the nasal cannula. Move your hand in front of the prongs. If you cannot feel the oxygen flow, check the cannula connections for leaks.

5.4 Alarm signals

The oxygen concentrator has the following alarm functions:

- 1) Pressure failure
- 2) Compressor failure
- 3) Low oxygen concentration
- 4) Overtemperature
- 5) Low oxygen flow rate
- 6) Power supply failure
- 7) Power failure
- 8) Start-up phase

5. OPERATION AND INSTALLATION

NOTE: All alarms on the device are low priority.

NOTE: All alarms are technical alarms.

When the oxygen concentrator is started, the blue, green and yellow indicators light up and the alarm sounds once to ensure that the alarm system is working properly. The blue and yellow indicators then go out. Five minutes after the oxygen concentrator is started, the oxygen sensor is working normally and controls the indicator lights according to the oxygen concentration value.

GB

Explanation of indicators and symbols

Symbol	Condition	Indicator lights	Alarm
OK	System is in good condition: Oxygen concentration $\geq$ 82%	green	---
	1) Oxygen concentration < minimum concentration (warm-up phase) 2) Oxygen concentration < 82%	yellow	Alarm
	System error (pressure error; compressor failure; overtemperature; low oxygen flow rate)	yellow	Alarm
	Power failure; power failure alarm	yellow	Alarm
	Audio alarm stopped	blue	---

5. OPERATION AND INSTALLATION

Description of alarm conditions

1. The oxygen concentration is below the minimum concentration during the start-up phase. The yellow lamp lights up and the total time is shown on the display. The device is in the warm-up phase. Wait 3 minutes. If the alarm persists, notify the supplier immediately.

2. The oxygen concentration is higher than 82%. The green indicator light is on and the total time is shown on the display. Normal operation.

3. Oxygen concentration below 82%. The yellow light illuminates, the alarm sounds and the control panel displays the total time. Contact your supplier immediately. You can continue to use the oxygen concentrator unless you receive different instructions from your supplier. Ensure that the backup oxygen is nearby.

NOTE: The oxygen concentrator reaches its most stable state after the warm-up phase (after approx. 15 minutes). The maximum and average delay of the low oxygen concentration alarm is 60 seconds.

4. In the event of a high/low pressure alarm, the yellow LED will illuminate, an alarm will sound, 'E1' or 'E2' will be displayed on the control panel, and the device will automatically switch off. If the device is still in operation, please switch it off manually. Use the reserve oxygen and contact the supplier immediately. The maximum and average delay of the 'E1' alarm system is less than 10 seconds. The maximum and average delay of the 'E2' alarm system is less than 5 seconds.

5. In the event of a compressor failure alarm, the yellow lamp will illuminate, an alarm will sound, 'E3' or 'E4' will appear on the control panel, and the unit will automatically shut down. If the unit is still in operation, please switch it off manually, use reserve oxygen, and contact the supplier immediately. The maximum and average delay of the alarm system due to a compressor failure is less than 10 seconds.

6. In the event of an overtemperature alarm, the yellow lamp will illuminate, an alarm will sound, 'E5' will appear on the control panel, and the device will automatically switch off. Please disconnect the power supply immediately, use reserve oxygen, and contact the supplier without delay. The maximum and average delay of the overtemperature alarm system is less than 10 seconds.

7. In the event of a low oxygen flow alarm, the yellow LED will illuminate, an alarm will sound, the word 'LL' will appear on the control panel, and the device will automatically switch off. Please disconnect the power supply immediately, use reserve oxygen, and contact your supplier immediately. The maximum and average delay of the low oxygen flow alarm is 32 seconds.

8. When a power failure alarm is triggered, the yellow LED lights up, the alarm sounds and the word 'E7' appears on the control panel. Please check the power supply. The maximum and average delay of the power failure alarm system is less than 10 seconds.

9) In the event of a power failure, the yellow LED lights up, the alarm sounds, nothing is displayed and the device switches off. Please check the power supply.

5. OPERATION AND INSTALLATION

Audio Alarm Pause Function: If the oxygen concentrator triggers an alarm, press the Audio Alarm Pause button to temporarily turn off the audible alarm. The blue light will illuminate. After 2 minutes or when pressed again, the alarm will resume and the blue light will turn off. The alarm sound is paused for 2 minutes and the oxygen concentrator returns to the alarm state after 2 minutes.

Operator position: Maximum distance of 1 metre from the oxygen concentrator.

GB

Alarm limits

Alarm	Alarm limits
High pressure	The pressure is more than 240 kPa.
Low pressure	The pressure is less than 20 kPa.
Compressor High current	The current is greater than 4.0 A(AC)
Compressor Low current	The current is equal to 0 A(AC)
Excessive temperature	The temperature of the air around the sensor is greater than 68°C.
Low oxygen concentration	The oxygen concentration is less than 82%.
Low flow rate	The flow rate is less than 0.3 l/min.
Power failure	The voltage is less than 185 ± 5 V (AC)
Voltage drop	The voltage is 0 V (AC).

NOTE: The alarm value is detected by the sensor.

5. OPERATION AND INSTALLATION

5.5 Timer function

This oxygen concentrator has a timer function that can be set by the user in the range of 0 to 6 hours.

When it starts working, the display shows „TIMING ----“, which means that the timer function is closed. The timer continues to run until the power supply is interrupted.

Press the „⏻“ button once to extend the timer by 30 minutes. If the button is held down for longer than 1,5 seconds, the timer will increase continuously. The oxygen concentrator switches off automatically and the display shows „TIMING ----“ and the total hours when the time measurement is complete. Reset the time measurement function.

5.6 Symbols

Symbol	Description	Symbol	Description
	Alternating current		Caution/Warning
	Class II device		Type BF application part
	OFF (disconnected from the mains)		ON (connected to the mains)
	Stack limit		Operation and storage only upright
	Temperature limit		Moisture limitation
	No smoking		No open flames: Fire, open ignition sources and smoking prohibited
	Keep dry		Fragile
	See operating instructions		Manufacturer

5. OPERATION AND INSTALLATION

Symbol	Description	Symbol	Description
	Atmospheric pressure limitation		European representative
	Audio alarm stopped		Alarm
	Standby		Time setting
	Brightness		Date of manufacture
	serial number		Medical device
	MR Unsure This device poses a significant safety risk to patients and staff if it is brought into the MR (magnetic resonance) environment!		
IP21	unsafe The first digit '2': protection against solid foreign objects with a diameter ≥ 12.5 mm $\varnothing$. The second digit '1': protection against vertically falling water drops.		

GB

5.7 Adjusting the brightness

Press the button „☀“, adjust the brightness of the display panel.

5.8 Switch off

While using the device, the user can press the „⏻“ button to stop/start the oxygen supply. First, remove the nasal cannula from the oxygen outlet, press the power switch to the „O“ position to turn off the oxygen concentrator, and then disconnect the power supply.

5. OPERATION AND INSTALLATION

5.9 Accessories

- The maximum oxygen flow rate for the accessories must not exceed 10 l/min. The maximum pressure for the accessories must not exceed 150 kPa.
 - This oxygen concentrator, its parts and accessories are designed for use at specific flow rates.
 - Incompatible parts or accessories may affect performance.
 - The responsible organisation is responsible for checking the compatibility of the oxygen concentrator and all parts or accessories used for oxygen therapy before use with the patient.
- ⚠ WARNING: Before and during oxygen therapy, only use water-based lotions or ointments that are compatible with oxygen.
- ⚠ CAUTION: Never use petroleum- or oil-based lotions or ointments to avoid the risk of fires and burns.
- ⚠ CAUTION: Correct placement and positioning of the nasal cannula in the nose is crucial for the amount of oxygen delivered to the patient's respiratory system.
- ⚠ CAUTION: The nasal cannula is intended for single use only and should be used immediately after opening the packaging and then disposed of.
- ⚠ CAUTION: The nasal cannula must not be used if the packaging has been damaged prior to use. Reusing the cannula may increase the risk of reinfection.
- ⚠ CAUTION: Using a nasal cannula other than the recommended one, e.g. a paediatric cannula for adult patients, may impair the effectiveness of oxygen therapy.
- ⚠ CAUTION: Only use accessories supplied or recommended by Yuwell to ensure compatibility between the device and accessories.
- ⚠ CAUTION: Recommended nasal cannula: PVC for adults, 2 m long, manufactured by JI-ANGSU WEIKANG JIEJING MEDICAL APPARATUS CO.
- ⚠ CAUTION: The fire-safe valve is a thermal fuse that extinguishes a fire in an oxygen tube and stops the flow of oxygen if the tube is accidentally opened. This user manual contains information about your concentrator and serves as a reference guide for using the concentrator. Please read it carefully and keep it in a safe place.

6. MAINTENANCE AND CARE

⚠ CAUTION: Before servicing the oxygen concentrator, it must first be disconnected from the mains power supply to avoid electric shock.

⚠ CAUTION: Under normal conditions and in the event of a single fault, the housing, humidifier and nasal cannula may become contaminated with bodily fluids or exhaled gases. To reduce the risk of infection, please perform maintenance on a routine basis.

NOTE: In locations with high dust or soot levels, maintenance may need to be performed more frequently.

NOTE: After cleaning and disinfecting the oxygen concentrator, parts or accessories, please wrap them in plastic bags and store them in a dry environment until the next use.

GB

6.1 Cleaning

⚠ CAUTION: Do not disassemble the outer casing of the oxygen concentrator.

⚠ WARNING: Liquids damage the internal components of the oxygen concentrator and its equipment. To avoid damage or injury from electric shock:

- Switch off the concentrator and disconnect the power cord before cleaning.
- Ensure that no cleaning agents enter the air inlet and outlet openings.
- Do NOT spray or apply cleaning agents directly onto the housing.
- Do NOT spray the product with a water hose.
- Do NOT immerse the device or accessories in liquid.

Clean the outer casing once a month as follows:

- 1) Use a cloth or sponge with a mild detergent or warm soapy water to clean the outer casing.
- 2) Allow the oxygen concentrator to air dry or use a dry towel before putting the oxygen concentrator into operation.

⚠ CAUTION: After cleaning and before drying, the oxygen concentrator must be disinfected as follows before it is made available to a new patient:

Wipe the outer casing with a cloth or sponge soaked in 70% to 80% medical alcohol.

6. MAINTENANCE AND CARE


6.2 Clean or replace filters

Please clean or replace the filters in good time, as this is very important for protecting the compressor and extending the service life of the oxygen concentrator.

Disassemble the filter:

Remove the filter cover to remove the filter screen. (Figure 9)

Fig. 9: Removing the filter

Clean the filter

- 1) Clean the filter with a mild detergent or warm soapy water and then rinse it thoroughly with clean water.
- 2) Dry the filter thoroughly before reinstalling it.
- 3) The filter must be cleaned or replaced once a month or as needed.

⚠ CAUTION: Do not operate the oxygen generator without a filter installed or if the filter is wet. Doing so may cause permanent damage to the oxygen concentrator.

6.3 Clean humidifiers


Disassemble the humidifier

Turn the humidifier bottle counterclockwise to open the humidifier and remove the humidifier hose and filter element. (Figure 10, Figure 11)

Abb. 10

Abb. 11

6. MAINTENANCE AND CARE

Clean the humidifier

Clean the humidifier monthly as follows to reduce limescale deposits and eliminate possible bacterial contamination:

- 1) Clean the parts of the humidifier with a mild detergent or warm soapy water and then rinse them thoroughly with clean water.
- 2) Allow to air dry thoroughly.

⚠ CAUTION: To limit bacterial growth, allow the humidifier to air dry thoroughly after cleaning when not in use. Fill the humidifier with clean water daily before use.

⚠ CAUTION: After cleaning and before drying, the oxygen concentrator must be disinfected as follows before it is made available to a new patient: Place the parts of the humidifier in 70% to 80% medical alcohol and leave them covered to soak for 30 minutes for disinfection.

⚠ The humidifier chamber is made of PP (polypropylene). It is suitable for prolonged disinfection in water at approximately 100 °C (boiling point). However, do not expose it to temperatures above 120 °C without water, as this may cause damage.

6.4 Cleaning the fire safety valve

Clean the fire safety valve monthly as follows:

- 1) Clean the Firesafe valve with a mild detergent or warm soapy water and then rinse thoroughly with clean water.
- 2) Allow it to dry thoroughly in the air.


⚠ CAUTION: After cleaning and before drying, the oxygen concentrator must be disinfected as follows before it is made available to a new patient: Immerse the fire safety valve in 70% to 80% medical alcohol and leave it covered to soak for 30 minutes for disinfection.

6.5 Replace fuse tube

Switch off the appliance, loosen the screws and separate the front and rear casing. Unscrew the fuse box anticlockwise and replace the fuse tube. (Figure 12)

Fuse type: T5AH250V, $\Phi 5 \times 20$

Fig. 12

6. MAINTENANCE AND CARE

6.6 Check the alarm system

The alarm system must be checked regularly: After the oxygen concentrator has been in operation for 5 minutes, set the flow meter to less than 0.3 l/min. After approx. 30 seconds, the yellow LED lights up, the alarm sounds, the word 'LL' appears on the control panel and the device switches off. Press the 'Audio Alarm Stop' button to turn off the alarm; the blue light will illuminate. Pressing the 'Audio Alarm Mute' button again will sound the alarm again and the blue light will turn off.

The methods for checking the alarm system's function for each of the alarm states are listed in the Technical Manual (document no.: 161849).

6.7 Instructions for preparation and reprocessing

- To prevent injury from infection or damage to the oxygen concentrator, only qualified personnel should clean and disinfect the oxygen concentrator and its accessories for use with multiple patients.
- Follow the instructions below to prevent possible infection with pathogens between patients due to contamination of components or accessories. If necessary, preventive maintenance should also be performed at this time.

- 1) Replace the nasal cannula, humidifier and oxygen tubing.
- 2) Check whether the oxygen concentrator is damaged or needs to be repaired.
- 3) Perform all procedures in the 'Maintenance' section.
- 4) Ensure that the oxygen concentrator is functioning normally and that all alarms are in good working order.
- 5) Before handing over to a new patient, ensure that the delivery includes the oxygen concentrator and this manual.

 **CAUTION:** Only the device itself may be reused. Humidifiers, oxygen tubing and nasal cannulas must always be replaced with new components.

7. TROUBLESHOOTING

Use the following table to take action if the oxygen concentrator indicates an unusual condition.

Problem	Potential Cause	Solution
The oxygen concentrator is not functioning, the yellow light is illuminated, the alarm is sounding, no display.	1) Poor contact between the mains cable plug and the socket.	1) Plug the power cord firmly into the power outlet.
	2) The socket has no power output.	2) Switch to a power outlet with power output.
	3) Insufficient power supply at the socket.	3) Do not use extension cables. Connect the oxygen concentrator to another power outlet.
	4) Fuse defective.	4) Replace the safety tube.
	5) If the oxygen concentrator still does not work, please contact the supplier.	
The oxygen concentrator is working and the operating noise is normal, the flow rate can be adjusted, but there is no oxygen supply or only a weak supply. Check the oxygen supply hose for kinks or blockages.	1) Air leak between humidifier bottle and cap.	1) Replace the humidifier bottle and cap and tighten them securely.
	2) The safety valve of the humidifier is open.	2) Shake the humidifier gently to close the safety valve.
	3) Air leak between humidifier and oxygen outlet..	3) Reinsert the humidifier.
	4) The accessories (nasal cannula, mask, humidifier, etc.) are leaking.	4) Replace the leaking accessory.
	5) If the problem persists, please contact the supplier.	

GB

7. TROUBLESHOOTING

Problem	Probable cause	Solution
The oxygen concentrator is working, but the yellow light is on and the alarm is sounding.	1) Oxygen concentration < 82%	1) Clean or replace the filter.
	2) Oxygen flow rate above the maximum recommended flow rate: 5 L/min.	2) Adjust the flow rate only on the advice of your doctor.
	3) If the problem persists, the device can still be used, but please contact the supplier in advance.	
The oxygen concentrator is not working, the yellow light is on, the alarm is sounding, and the word 'E1' is displayed on the control panel.	1) The system pressure is too low.	1) Clean or replace the filter.
	2) If the problem persists, stop using the device and contact the supplier immediately.	
The oxygen concentrator is not functioning; the yellow light illuminates, the alarm sounds, and the word 'E2' appears on the control panel.	1) The system pressure is too high.	1) Stop using the device and contact the supplier immediately.
The oxygen concentrator is not functioning; the yellow light illuminates, the alarm sounds, and the word 'E3' appears on the control panel.	1) The compressor circuit is open.	1) Stop using the device and contact the supplier immediately.

7. TROUBLESHOOTING

Problem	Probable cause	Solution
The oxygen concentrator is not functioning, the yellow light is illuminated, the alarm is sounding, and the word 'E4' is displayed on the control panel.	1) The compressor circuit is short-circuited.	1) Stop using the device and contact the supplier immediately.
The oxygen concentrator is not functioning, the yellow light is illuminated, the alarm is sounding, and the word 'E5' is displayed on the control panel.	1) The temperature inside the oxygen concentrator is too high.	1) Stop using the device and contact the supplier immediately.
The oxygen concentrator is working, but the yellow light is on, the alarm is sounding, and the word 'E7' is displayed on the control panel.	1) The voltage is too low.	1) Change the power source to achieve normal voltage.
	2) Stop using the device and contact the supplier immediately.	
The oxygen concentrator is not functioning, the yellow light is illuminated, the alarm is sounding, and the word 'LL' is displayed on the control panel.	1) The oxygen flow rate is too low.	1) Turn the flow meter knob counterclockwise to increase the flow.
	2) If the problem persists, stop using the device and contact the supplier immediately.	


GB

⚠ CAUTION: If you encounter any other problems, please switch the concentrator OFF first, use your oxygen reserve and contact the supplier immediately.


⚠ CAUTION: The steps for disassembling the oxygen concentrator (for service personnel only) are described in the Technical Manual (document no.: 161849).

8. FURTHER IMPORTANT INFORMATION

8.1 Functional diagram of oxygen production

8.2 Electrical function diagram

8. FURTHER IMPORTANT INFORMATION

8.3 Package contents

- 1x Oxygen concentrator
- 1x Manual
- 1x Filter
- 2x Safety tube: T5AH250V, $\Phi 5 \times 20$

Information for accessories

Name	Manufacturer	Type	Technical specifications
Fire safety valve	JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD	YY-ZYJTY-10-00	ABS, outer diameter of the plug: 7 mm

8.4 Disposal of equipment

The person responsible must contact the local authorities to determine the proper disposal method for the oxygen concentrator and accessories. Special attention must be paid to the disposal of molecular sieves.

8.5 Language

This operating manual is written in German.

9. EMC INFORMATION

- ⚠ • WARNING: Use the device well away from HF SURGICAL EQUIPMENT and the HF-shielded room of an MRI SYSTEM in hospitals, where the intensity of EM INTERFERENCE is high.
- ⚠ • WARNING: Avoid using this device next to or on top of other devices, as this may cause malfunctions. If such use is necessary, this device and the other devices should be observed to ensure that they are functioning normally.
- ⚠ • WARNING: The use of accessories, converters and cables not specified or supplied by the manufacturer of this device may result in increased electromagnetic emissions or reduced electromagnetic immunity of this device and may cause it to malfunction.
- ⚠ • WARNING: Portable RF communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) to any part of the 9F-5B oxygen concentrator, including the cables specified by the manufacturer. Failure to do so may affect the performance of the device.
- ESSENTIAL PERFORMANCE: The oxygen concentration in the delivered gas, both under normal conditions and in the event of a single fault, within the performance levels specified in the instructions for use, or the generation of an alarm condition: power failure, technical alarm condition, low oxygen concentration, technical alarm condition, malfunction, technical alarm condition, technical alarm condition during the start-up phase.
- ⚠ • CAUTION: If the oxygen concentrator is not functioning normally or an alarm condition occurs, the user should attempt to move the oxygen concentrator to another area to determine if the problem is due to electromagnetic interference from other devices in the vicinity.

Table 1:
Manufacturer's instructions and explanation – electromagnetic immunity

Problem	Basic EMC standard or test method	Interference immunity level
Electrostatic discharge	IEC 61000-4-2	±8 kV Contact ±2 kV, ±4 kV, ±8 kV, ±15 kV Air
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz - 2,7 GHz 80 % AM bei 1 kHz
Nominal power Frequency Magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz
Fast electrical transients/bursts	IEC 61000-4-4	±2 kV 100 kHz Wiederholungsfrequenz

9. EMC INFORMATION

Problem	Basic EMC standard or test method	Interference immunity level
Overvoltages Line-to-line	IEC 61000-4-5	±0,5 kV, ±1 kV
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V/m 0,15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Voltage dips	IEC 61000-4-11	0 % U ; 0,5 Cycle T Bei 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0 % U ; 1 Cycle T and 70 % U; 25/30 Cycles T Single-phase: at 0°
Power outages	IEC 61000-4-11	0 % U; 250/300 Cycles

GB

Table 2:
Test specifications for the IMMUNITY OF ENCLOSURE CONNECTORS for wireless RF communication devices

Test frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 to 470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	28
710	704 to 787	LTE Band 13,17	Pulse modulation 217 Hz	9
745				
780				


9. EMC INFORMATION

Test frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITÄT TESTLEVEL (V/m)
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				
NOTE: If necessary to achieve IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME DEVICE or ME SYSTEM may be reduced to 1 metre. The test distance of 1 metre is permitted in accordance with IEC 61000-4-3.				

Table 3:
Instructions and manufacturer's declaration – electromagnetic radiation

Problem	Conformity
Conducted and radiated RF EMISSIONS CISPR 11	Group II, Class B
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Conform

GB

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD.
No.1 Baisheng Road Development Zone, Danyang, Jiangsu 212300 CHINA
TEL: 86-511-86900833
www.yuwell.com

Metrax GmbH
Rheinwaldstr. 22, 78628 Rottweil, GERMANY

yuwell  0123

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